

ANNEX I: Common list of rapid antigen tests¹

As agreed by Member States on [...]

Disclaimer: *This list was agreed by the HSC based on a proposal by the Technical Working Group on COVID-19 Diagnostic Tests. Experts participating in the Technical Working Group strongly recommend that use of rapid antigen tests is primarily intended for preliminary testing for SARS-CoV-2 infection in symptomatic patients, and note that rapid antigen tests should in particular be used in the specific contexts and circumstances referred to by the Commission Recommendation (EU) 2020/1743 of 18 November 2020 and the technical guidance by ECDC on 19 November 2020. The content of the common list is based on the clinical performance data and information that is available at this moment in time. Updates to the common list are based on the criteria as described in Council Recommendation 2021/C 24/01 as well as the additional criteria and definitions agreed by the Technical Working Group on 14 September 2021 (tbc). The Medical Device Coordination Group Guidance on performance evaluation of SARS-CoV-2 in vitro diagnostic medical devices², envisaged to form the basis for common specifications to be adopted according to Article 9 of Regulation (EU) 2017/746, has been taken into consideration in this regard.*

Questions to TWG:

- Please review the entries below. In particular, for those tests where the study design is listed as ‘unknown’, please assess these cases, collect more information and evaluate if and how these evaluation studies and their results should be listed.
- In case any new data or results are known from evaluation studies (matching the criteria and definition) not yet listed, please provide these details.
- Which tests included in the list are ‘twin tests’?

¹ This is the list of rapid antigen tests as referred to in Article 3 of the Regulation (EU) 2021/953 of the European Parliament and of the Council of 14 June 2021 on a framework for the issuance, verification and acceptance of interoperable COVID-19 vaccination, test and recovery certificates (EU Digital COVID Certificate) to facilitate free movement during the COVID-19 pandemic, OJ L 211, 15.6.2021, p. 1–22.

² https://ec.europa.eu/health/sites/default/files/md_sector/docs/mdcg_2021-21_en.pdf

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance Data by manufacturer	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
AAZ-LMB	COVID-VIRO® Rapid antigen test COVID-19	1833	<i>Study design unknown</i>						
			BE: 96.6% sensitivity, 100% specificity, NP swab FR: >95% sensitivity, 100% specificity SI: 96.6% sensitivity, 100% specificity, NP swab	96.6% sensitivity 100% specificity	BE, FR, SI	CH	FR CH		10 May 2021
Abbott Rapid Diagnostics	Panbio™ COVID-19 Ag Rapid Test	1232	<i>Prospective clinical field study</i>						
			FIND evaluation studies DE (10 Dec 2020) 1108 samples, NP swab Clinical sensitivities: - Days ≤ 7: 90.8%; - Ct ≤ 33: 88.3%; - Ct ≤ 25: 95.8%; Clinical specificity: 99.9%						
			CH (10 Dec 2020) 535 samples, NP swab Clinical sensitivities: - Days ≤ 7: 85.6%; - Ct ≤ 33: 89.7%; - Ct ≤ 25: 96.8%; Clinical specificity: 100%	91.4% sensitivity 99.8% specificity NP swab (Ct ≤ 33)	AT, BE, BG, CY, CZ, DE ^[2] , DK, EE, EL, ES, FR ^[1] , HR, IT, LT, LV, MT, NL ^[5] , PL, PT, RO, SE, SK	CH, ME, MK, NO, UK, UA	DE ^[2] , ES, FI, NL ^[5] , PT CH, NO	CY, ES, HR, HU, IE, LU, SE	17 February 2021
			India (25 June 2021) 526 samples, NP swab Clinical sensitivities: - Days ≤ 7: 61.3%-100%; - Ct ≤ 33: 74.2%-86.7%; - Ct ≤ 25: 91.9%-100%; Clinical specificity: 100%	98.1% sensitivity 99.8% specificity Nasal swab (Ct ≤ 33)					
			<i>Retrospective in vitro studies</i>						
			FI: Validated in several laboratories (studies not published), meeting criteria.						
			<i>Study design unknown</i>						
			DE: 91.4% sensitivity 99.8% specificity, NP swab; 98.1% sensitivity,						

³ As registered in and used by the JRC database, see: <https://covid-19-diagnostics.jrc.ec.europa.eu/>.

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance Data by manufacturer	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
			99,8 specificity, Nasal swab BE ^[6] : Small-scale head-to-head comparison of 5 RATs in Belgian hospital lab. Panbio overall sensitivity (Ct range 14,6 – 35,5): 45/57 samples (79%). Sensitivity for Ct≤25: 17/18 samples. Overall specificity 100%.						
Acon Biotech (Hangzhou) Co., Ltd	SARS-CoV-2 Antigen Rapid Test	1457	<i>Retrospective in vitro studies</i>	96.9% sensitivity Nasal swab	DE ^[2] , FR, PT		DE ^[2]		14 July 2021
			DE : Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 94,1% at <Ct25) + Manufacturer specificity: 99.54%						
ACON Laboratories, Inc.	Flowflex SARS-CoV-2 Antigen Rapid Test	1468	<i>Prospective clinical field study</i>	96.9% sensitivity Nasal swab	AT, BE, DE ^[2] , LT, LV, SI		DE ^[2]		10 May 2021
			FIND evaluation CH (9 June 2021) 279 samples, nasal swab Clinical sensitivities: - Days ≤ 7: 92.2%; - Ct ≤ 33: 98.3%; - Ct ≤ 25: 100%; Clinical specificity: 99.5%						
			<i>Retrospective in vitro studies</i>						
			DE : Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 94,1% at <Ct25) + Manufacturer specificity: 98,7%						
			<i>Study design unknown</i>						
			BE : 96.9% sensitivity, 99.5% specificity, NP swab						
AESKU.DIAGNOSTICS GmbH & Co, KG	AESKU.RAPID SARS-CoV-2	2108	<i>Study design unknown</i>	96% sensitivity 98% specificity NP swab	AT, DE ^[2] , SI		DE ^[2]		10 May 2021
			DE : 96% sensitivity, 98% specificity SI : 96% sensitivity, 98% specificity, Nasal swab						
Affimedix Inc.	TestNOW® - COVID-19 Antigen Test	2130	<i>Retrospective in vitro studies</i>	96.1% sensitivity 99.4% specificity NP swab	DE ^[2]		DE ^[2]		10 May 2021
			DE : Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 99,4%						
AMEDA	AMP Rapid Test SARS-	1304	<i>Retrospective in vitro studies</i>	97.3% sensitivity	AT, BG,	CH, UA	DE ^[2]	HR	17 February

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance Data by manufacturer	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
Labordiagnostik GmbH	CoV-2 Ag		DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 100% <i>Study design unknown</i> BE: 97.3% sensitivity, 100% specificity, NP swab SI: 97.3% sensitivity, 100% specificity, NP swab	NP swab 97.3% sensitivity Nasal swab 100% specificity	DE ^[2] HR, SI		CH		2021
Anbio (Xiamen) Biotechnology Co., Ltd	Rapid COVID-19 Antigen-Test (colloidal Gold)	1822	<i>Study design unknown</i> DE: 99.27% sensitivity, 100% specificity		AT, DE ^[2]		DE ^[2]		10 May 2021
Anhui Deep Blue Medical Technology Co., Ltd	COVID-19 (SARS-CoV-2) Antigen Test Kit (Colloidal Gold)	1736	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: >99% <i>Study design unknown</i> BE: 95% sensitivity, 99% specificity, NP/OP swab	Nasal swab: 96,4% sensitivity, 99,8% specificity NP swab: 95,7% sensitivity, 99,3% specificity OP swab: 96,4% sensitivity, 99,8% specificity	BE, DE ^[2]	UK	DE ^[2]		10 May 2021
Anhui Deep Blue Medical Technology Co., Ltd	COVID-19 (SARS-CoV-2) Antigen Test Kit (Colloidal Gold) – Nasal swab	1815	<i>Study design unknown</i> DE: 96,4 % sensitivity, 99,8 % specificity	96.4 % sensitivity 99.8 % specificity Nasal swab	DE ^[2]		DE ^[2]		10 May 2021
ArcDia International Ltd	mariPOC SARS-CoV-2	768	<i>Study design unknown</i> FI: Meets the minimum performance requirements – see the report for details.	92% sensitivity 100% specificity	FI		FI		10 May 2021
ArcDia International Oy Ltd	mariPOC Respi+	2078	<i>Retrospective in vitro studies</i> FI: Validated in several laboratories (studies not published), meeting criteria.	100 % sensitivity 100 % specificity NP swab	FI, PT		FI		14 July 2021
ArcDia International Oy Ltd	mariPOC Quick Flu+	2079	<i>Retrospective in vitro studies</i> FI: Validated in several laboratories (studies not published), meeting criteria.	100 % sensitivity 100 % specificity NP swab	FI, PT		FI		14 July 2021
Artron Laboratories Inc.	Artron COVID-19 Antigen Test	1618	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 100%	96.67% (Nasal)sensitivity 91.67% (NP) sensitivity 100 % specificity Nasal/NP swab	DE ^[2]		DE ^[2]		14 July 2021
Asan Pharmaceutical	Asan Easy Test COVID-19	1654	<i>Study design unknown</i>		DE ^[2]		DE ^[2]		10 May 2021

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance Data by manufacturer	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
Co., Ltd	Ag		DE: 94.67% sensitivity, 97.71% specificity						
Assure Tech. (Hangzhou) Co., Ltd.	ECOTEST COVID-19 Antigen Rapid Test Device	770	<i>Retrospective in vitro studies</i>	92.5 % sensitivity 99.2 % specificity Nasal/NP/ OP swab	DE ^[2]		DE ^[2]		14 July 2021
			DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 95% at <Ct25) + Manufacturer specificity: 99.2%						
Assure Tech. (Hangzhou) Co., Ltd.	ECOTEST COVID-19 Antigen Rapid Test Device	2350	<i>Retrospective in vitro studies</i>	Sensitivity: 97.7%, Specificity: 99.1% NP and OP swab	CZ, DE ^[2]		DE ^[2]		23 July 2021
			DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 95% at <Ct25) + Manufacturer specificity: 99.1%						
Atlas Link Technology Co. Ltd.	NOVA Test ® SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold Immunochromatography)	2010	<i>Study design unknown</i>	98.5 % sensitivity 99.4 % specificity Nasal/OP swab	AT, DE ^[2] , SI	CH	DE ^[2] CH		10 May 2021
			DE: 97.6% sensitivity, 99.2% specificity						
Avalun	Ksmart® SARS-COV2 Antigen Rapid Test	1800	<i>Retrospective in vitro studies</i>	Clinical Sensitivity: 93.18 % Clinical Specificity: 99.32 % NP swab	DE ^[2]		DE ^[2]		7 July 2021
			DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 94,1% at <Ct25) + Manufacturer specificity: 99,32%						
AXIOM Gesellschaft für Diagnostica und Biochemica mbH	COVID-19 Antigen Rapid Test	2101	<i>Study design unknown</i>	98% sensitivity 100% specificity NP/Nasal swab	DE ^[2]		DE ^[2]		10 May 2021
			DE: 98.1% sensitivity, 100% specificity						
Azure Biotech, Inc.	COVID-19 Antigen Rapid Test Device	1906	<i>Study design unknown</i>	95% sensitivity 99.2% specificity NP swab	DE ^[2]		DE ^[2]		10 May 2021
			DE: 94.3% sensitivity, 99.1% specificity						
Becton Dickinson	BD Veritor™ System for Rapid Detection of SARS CoV 2	1065	<i>Prospective clinical field study</i>	Clinical Sensitivity: 91.1 % Clinical Specificity: 99.6 % Nasal swab	NL		NL		7 July 2021
			NL: Independent field study in symptomatic individuals - sampling was Nasal mid-turbinate and OP swab. Sensitivity overall: 79.5% - Sensitivity Ct<30: 93.2% - Specificity overall: 99.8%						
Beijing Hotgen	Novel Coronavirus 2019-	1870	<i>Retrospective in vitro studies</i>	97.1% sensitivity	AT, BE,		DE ^[2]		10 May 2021

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance Data by manufacturer	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
Biotech Co., Ltd	nCoV Antigen Test (Colloidal Gold)		DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 99.76% <i>Study design unknown</i> BE: 98.6% sensitivity, 100% specificity, NP Swab 97.3% sensitivity, 99.2% specificity. OP swab SI: 96.6% sensitivity, 99.8% specificity, NP swab	99.76% specificity	DE ^[2] , RO, SI				
Beijing Jinwofu Bioengineering Technology Co.,Ltd.	Novel Coronavirus (SARS-CoV-2) Antigen Rapid Test Kit	2072	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 100%	96.88 % sensitivity 100 % specificity Nasal/ NP/ OP swab	DE ^[2]		DE ^[2]		14 July 2021
Beijing Lepu Medical Technology Co., Ltd	SARS-CoV-2 Antigen Rapid Test Kit	1331	<i>Study design unknown</i> BE: 92% sensitivity, 99.3% specificity, Nasal DE: 92.0% sensitivity, 99.26% specificity SI: 92% sensitivity, 99.2% specificity, NP	92% sensitivity Nasal swab	AT, BE, DE ^[2] , SI, RO	UA	DE ^[2]		17 February 2021
Beijing Wantai Biological Pharmacy Enterprise Co., Ltd	Wantai SARS-CoV-2 Ag Rapid Test (FIA)	1484	<i>Study design unknown</i> DE: 96.6% sensitivity, 96.9% specificity	96.6% sensitivity, Nasal swab	DE ^[2]		DE ^[2]		17 February 2021
Beijing Wantai Biological Pharmacy Enterprise Co., Ltd	Wantai SARS-CoV-2 Ag Rapid Test (colloidal gold)	1485	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 99%	96.1 % sensitivity 99% specificity Nasal swab	DE ^[2]		DE ^[2]		14 July 2021
BioGnost Ltd	CoviGnost AG Test Device 1x20	2247	<i>Retrospective in vitro studies</i> HR: 300 NP samples (retrospective), symptomatic (<7 dps): 200 PCR+ samples (range Ct 16-30), Ct<30: sensitivity 96.5% 100 PCR- samples: specificity 100%	Sensitivity: 96%, Specificity: 99% NP swab	HR		HR		23 July 2021
BIOHIT HealthCcare (Hefei) Co., Ltd.	SARS-CoV-2 Antigen Rapid Test Kit (Fluorescence Immunochromatography)	Yes (1286)	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 98.9%	Sensitivity: 96.77% Specificity: 98.9% NP/OP swab	DE ^[2]		DE ^[2]		23 July 2021
BioMaxima SA	SARS-CoV-2 Ag Rapid Test	Yes	<i>Study design unknown</i>	Sensitivity: 95%	PL		PL		23 July 2021

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance Data by manufacturer	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
		(2035)	PL: Diagnostic sensitivity: 93.43% (95% CI: 91.61%~97.19%); diagnostic specificity: 97.75%, manufacturer specificity: 99.1%	Specificity: 99% NP Swab					
Biomerica Inc.	Biomerica COVID-19 Antigen Rapid Test (nasopharyngeal swab)	1599	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 99,7%	Clinical Sensitivity: 94.7 % Nasal/NP swab	DE ^[2]		DE ^[2]		7 July 2021
BIONOTE	NowCheck COVID-19 Ag Test	1242	<i>Prospective clinical field study</i> FIND evaluation Brazil (20 April 2021) 400 samples, NP swab Clinical sensitivities: - Days ≤ 7: 92.2%; - Ct ≤ 33: 91.4%; - Ct ≤ 25: 94.8%; Clinical specificity: 97.3% Brazil (30 March 2021) 218 samples, Nasal/NP swab. Clinical sensitivities: - Days ≤ 7: 92.5% (N/NP); - Ct ≤ 33: 97.2% (N/NP); - Ct ≤ 25: 100% (N/NP); Clinical specificity: 98.6% <i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 98,6%	Clinical Sensitivity: 90.91 % Clinical Specificity: 99.43 %	DE ^[2]		DE ^[2]		7 July 2021
BIO-RAD	CORONAVIRUS AG RAPID TEST CASSETTE	2031	<i>Study design unknown</i> ES: NP swab: sensitivity 98,3%; specificity 99,6% (119 positive samples, 746 negative samples) Nasal swab: sensitivity 97,2%; specificity 100% (109 positive samples, 128 negative samples)	Clinical Sensitivity: 98 % (NP: 98,32% / Nasal: 97,25%) Clinical Specificity: 99 % (NP: 99,6% / Nasal: 100%)	ES		ES		7 July 2021
BIOSYNEX S.A.	BIOSYNEX COVID-19 Ag BSS	1223	<i>Prospective clinical field study</i> NL: Independent field study, mainly symptomatic individuals, sensitivity Ct≤30: 96.0%; specificity overall: 100% <i>Retrospective in vitro studies</i>	96% sensitivity, 100% specificity, NP swab	AT, BE, DE ^[2] , DK, FR, NL ^[5] , PT	CH	DE ^[2] , NL ^[5] , CH	DK	17 February 2021

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance Data by manufacturer	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
			<p>DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 100%</p> <p><i>Study design unknown</i></p> <p>BE^[6]: Small-scale head-to-head comparison of 5 RATs in Belgian hospital lab. Biosynex overall sensitivity (Ct range 14,6 – 35,5): 52/58 samples (89,7%). Sensitivity for Ct≤25: 18/18 samples. Overall specificity only 46,2%, but this is probably linked to the use of transport medium instead of the swab included in the kit.</p>						
BIOSYNEX SA	BIOSYNEX COVID-19 Ag+ BSS	1494	<p><i>Prospective clinical field study</i></p> <p>FR: Validation study data: 125 positive and 118 negative samples; sensitivity 96%, specificity: 99%</p>	Clinical Sensitivity: 97.5 %	FR		FR		7 July 2021
BIOTEKE CORPORATION (WUXI) CO., LTD	SARS-CoV-2 Antigen Test Kit (colloidal gold method)	2067	<p><i>Retrospective in vitro studies</i></p> <p>DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 95% at <Ct25) + Manufacturer specificity: 99.28%</p>	96.49 % sensitivity 99.28 % specificity OP/NP swab	DE ^[2]		DE ^[2]		14 July 2021
Biotical Health S.L.U. BIOTICAL HEALTH S.L.U	biotical SARS-CoV-2 Ag Card	Yes (2013)	<p><i>Retrospective in vitro studies</i></p> <p>BE: Validation study 1: sensitivity 91.7% for Ct<25; Validation study 2: 94% for Ct<25. Manufacturer specificity: 99%</p>	Sensitivity: 96%, Specificity: 99% NP swab	BE		BE		23 July 2021
Boditech Med Inc	AFIAS COVID-19 Ag	Yes (1989)	<p><i>Prospective clinical field study</i></p> <p>NL: Independent field study in mild symptomatic (n= 427); overall sensitivity: 81.1% (106 PCR+), Ct <30: 96.4% (85 PCR+), PCR on NP+OP, Target antigen = nucleoprotein</p>	Sensitivity: 91.7%, Specificity: 98.7% NP swab	FR, NL		NL		23 July 2021
BTNX Inc	Rapid Response COVID-19 Antigen Rapid Test	1236	<p><i>Study design unknown</i></p> <p>DE: 94.55% sensitivity, 100% specificity</p>	90.2% sensitivity 100% specificity NP swab, NP swab, OP swab	AT, DE ^[2] , ES, SI		DE ^[2]		10 May 2021
CerTest Biotec	CerTest SARS-CoV-2 Card test	1173	<p><i>Prospective clinical field study</i></p> <p>ES: Ct ≤ 25, sensitivity: 94,0%; sensitivity for samples within the first 5 days after symptom onset: 84,8%</p>	92.9% sensitivity 99.6% specificity NP swab	ES, PT, SI		DE ^[2] , ES		17 February 2021
Core Technology	Coretests COVID-19 Ag	1919	<i>Study design unknown</i>	98.1% sensitivity	AT, DE ^[2] ,		DE ^[2]		10 May 2021

[illegible]

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance Data by manufacturer	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
			DE: Positive evaluation by Paul-Ehrlich-Institut (Sensitivity of 100% at <Ct25) + Manufacturer Specificity: 99,24%						
Eurobio Scientific	EBS SARS-CoV-2 Ag Rapid Test	1739	Prospective clinical field study	Clinical Sensitivity: 95.7 % Nasal swab	DE ^[2] , FR		DE ^[2] , FR		7 July 2021
			FR: Validation study data: 119 positive and 125 negative samples; sensitivity 93%, specificity: 99%						
			Retrospective in vitro studies						
			DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 94,1% at <Ct25) + Manufacturer specificity: 99,1%						
Fujirebio	ESPLINE SARS-CoV-2	2147	Prospective clinical field study	Clinical Sensitivity: 87.8 % (n=98, Ct<33) Clinical Specificity: 100 % NP swab	DE ^[2]		DE ^[2]		7 July 2021
			FIND evaluation DE (29 March 2021) 723 samples, NP swab Clinical sensitivities: - Days ≤ 7: 88.5%; - Ct ≤ 33: 87.8%; - Ct ≤ 25: 92.4%; Clinical specificity: 100%						
			Retrospective in vitro studies						
			DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 94,1% at <Ct25) + Manufacturer specificity: 99,13%						
GA Generic Assays GmbH	GA CoV-2 Antigen Rapid Test	Yes (1855)	Retrospective in vitro studies	Sensitivity: 97.059%, Specificity: 99.2% NP swab	DE ^[2]		DE ^[2]		23 July 2021
			DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 99.2%						
GenBody Inc	Genbody COVID-19 Ag Test	1244	Study design unknown	90% sensitivity 98% specificity NP/OP swab	DE ^[2]	UA	DE ^[2]		17 February 2021
			DE: 90% sensitivity 98% specificity						
Genrui Biotech Inc	SARS-CoV-2 Antigen Test Kit (Colloidal Gold)	2012	Retrospective in vitro studies	Sensitivity: 91.15% Specificity: 99.02% Nasal/NP/OP swab	DE ^[2]		DE ^[2]		7 July 2021
			DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 94,1% at <Ct25) + Manufacturer specificity: 99,02%						
GenSure Biotech Inc	GenSure COVID-19	1253	Retrospective in vitro studies	96.86% sensitivity,	DE ^[2]		DE ^[2]		10 May 2021

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance <i>Data by manufacturer</i>	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
	Antigen Rapid Test Kit		DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 94,1% at <Ct25) + Manufacturer specificity: 100%	100% specificity Nasal swab					
Getein Biotech, Inc	SARS-CoV-2 Antigen (Colloidal Gold)	1820	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 98.71%	97.06% sensitivity 98.71% specificity Nasal swab	AT, DE ^[2]		DE ^[2]		14 July 2021
Getein Biotech, Inc.	One Step Test for SARS-CoV-2 Antigen (Colloidal Gold)	2183	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 90% at <Ct30 and 100% at <Ct25)	97.06% sensitivity 98.71% specificity Nasal swab	DE ^[2]		DE ^[2]		16 June 2021
Goldsite Diagnostic Inc.	SARS-CoV-2 Antigen Kit (Colloidal Gold)	1197	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity 100% at <Ct25)		BE, BG, CY, FR, RO, SI, ES	UK	FR, DE ^[2] , ES		14 July 2021
Green Cross Medical Science Corp.	GENEDIA W COVID-19 Ag	1144	<i>Study design unknown</i> BE: 90.2% sensitivity, 100% specificity, NP swab DE: 90.1% sensitivity, 100% specificity	100% sensitivity 90.1% sensitivity NP swab, Anterior nasal swab	AT, BE, DE ^[2]		DE ^[2]		10 May 2021
Guangdong Hecin Scientific, Inc.	2019-nCoV Antigen Test Kit (colloidal gold method)	1747	<i>Study design unknown</i> DE: 96.6% sensitivity, 99.07% specificity	96.23% sensitivity Nasal swab	AT, DE ^[2]		DE ^[2]		10 May 2021
Guangdong Longsee Biomedical Co., Ltd.	COVID-2019-nCoV Ag Rapid TestDetection Kit(Immuno-Chromatography)	1216	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 99.5%	99.72% specificity NP/OP swab	DE ^[2]		DE ^[2]		14 July 2021
Guangdong Wesail Biotech Co. Ltd	COVID-19 Ag Test Kit	1360	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 98% <i>Study design unknown</i> SI: 90% sensitivity, 98% specificity, NP/Nasal swab	90% sensitivity 98% specificity Nasal swab	DE ^[2] , SI		DE ^[2]		17 February 2021
Guangzhou Decheng Biotechnology CO., Ltd	V-CHEK, 2019-nCoV Ag Rapid Test Kit (Immunochromatography)	1324	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 94,1% at <Ct25) + Manufacturer specificity: 99,5%	Clinical Sensitivity: 96.67 % Nasal swab	DE ^[2]		DE ^[2]		7 July 2021
Guangzhou Wondfo	Wondfo 2019-nCoV	1437	<i>Prospective clinical field study</i>		AT, BE, BG,	CH	DE ^[2]		10 May 2021

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance Data by manufacturer	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
Biotech Co., Ltd	Antigen Test (Lateral Flow Method)		FIND evaluation CH (25 Feb 2020) 328 samples, NP swab Clinical sensitivities: - Days \leq 7: 85.7%; - Ct \leq 33: 92.2%; - Ct \leq 25: 100%; Clinical specificity: 100% <i>Study design unknown</i> BE: 96.2% sensitivity, 99.7% specificity, NP/OP swab DE: 96.18 % sensitivity, 99.72% specificity		DE ^[2] , FR				
Hangzhou Lysun Biotechnology Co. Ltd	COVID-19 Antigen Rapid Test Device (Colloidal Gold)	2139	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 100%	96.46% sensitivity 100% specificity Nasal swab	DE ^[2]	CH	DE ^[2]		10 May 2021
Hangzhou AllTest Biotech Co., Ltd	COVID-19 Antigen Rapid Test	1257	<i>Study design unknown</i> DE: 93,40% sensitivity, 99,90% specificity	NP swab	AT, BE, BG, FR, SI, RO	CH	DE ^[2]	AT	10 May 2021
Hangzhou Clongene Biotech Co., Ltd	COVID-19 Antigen Rapid Test Casette	1610	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 94,4% at <Ct25) + Manufacturer specificity: 100%	Clinical Sensitivity: 91.4 % Clinical Specificity: 100 % NP swab	DE ^[2]		DE ^[2]		7 July 2021
Hangzhou Clongene Biotech Co., Ltd.	Covid-19 Antigen Rapid Test Kit	1363	<i>Study design unknown</i> BE: 91.4% sensitivity, 100% specificity, NP/OP swab DE: 91.4% sensitivity, 99.4% specificity SI: 91.4% sensitivity, 100% specificity, NP/OP swab	98.5% (Ct<33) sensitivity Nasal swab	AT, BE, DE ^[2] , FR, SI	CH	DE ^[2] CH	HR	17 February 2021
Hangzhou Clongene Biotech Co., Ltd.	COVID-19/Influenza A+B Antigen Combo Rapid Test	1365	<i>Study design unknown</i> DE: 97.7% sensitivity, 99.8% specificity	91% sensitivity 100% specificity NP swab	DE ^[2]		DE ^[2]		10 May 2021
Hangzhou Immuno Biotech Co., Ltd	Immunobio SARS-CoV-2 Antigen ANTERIOR NASAL Rapid Test Kit (minimal invasive)	1844	<i>Study design unknown</i> DE: 94.39% sensitivity 97.67% specificity	94% sensitivity 100% specificity Nasal swab, NP	DE ^[2]		DE ^[2]		10 May 2021
Hangzhou Immuno	SARS-CoV2 Antigen Rapid	2317	<i>Study design unknown</i>	Clinical Sensitivity	AT, DE ^[2]		DE ^[2]		10 May 2021

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance Data by manufacturer	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
Biotech Co., Ltd	Test		DE: 95.6% sensitivity, 100% specificity	98 % Clinical Specificity 100 %					
Hangzhou Laihe Biotech Co.	LYHER Novel Coronavirus (COVID-19) Antigen Test Kit (Colloidal Gold)	1215	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 94,1% at <Ct25) + Manufacturer specificity: 99,7%	Clinical Sensitivity: 95.07% % Clinical Specificity: 99.74% Nasal swab	AT, DE ^[2]	CH	DE ^[2]		10 May 2021
Hangzhou Testsea Biotechnology Co., Ltd.	Covid-19 Antigen Test Cassette	1392	<i>Study design unknown</i> DE: 97.6% sensitivity 98.4% specificity	92.1% sensitivity 98.1% specificity Nasal swab	DE ^[2]		DE ^[2]		10 May 2021
Healgen Scientific	Coronavirus Ag Rapid Test Cassette	1767	<i>Study design unknown</i> DE: 97.25% sensitivity, 100% specificity SI: 96.7% sensitivity, 99.2% specificity, NP/Nasal swab	80.6 % sensitivity 99.7% specificity NP swab	AT, DE ^[2] , NL ^[5] , SE, SI	CH	DE ^[2] , NL ^[5]	SE ^[3]	17 February 2021
Hubei Jinjian Biology Co., Ltd	SARS-CoV-2 Antigen Test Kit	Yes (1759)	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 99.3 %	Sensitivity: 98.02% Nasal Swab	DE ^[2]		DE ^[2]		23 July 2021
Humasis	Humasis COVID-19 Ag Test	1263	<i>Study design unknown</i> BE: 95.5% sensitivity, 100% specificity, NP swab DE: 95.5% sensitivity, 100% specificity SI: 95.5% sensitivity, 100% specificity, NP swab	95.3% sensitivity Nasal swab	AT, BE, BG, DE ^[2] , FR, HR, SE, SI		DE ^[2]	HR, SE	10 May 2021
Jiangsu Bioperfectus Technologies Co., Ltd.	Novel Corona Virus (SARS-CoV-2) Ag Rapid Test Kit	2107	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 99.15%	97.06 % sensitivity 99.15 % specificity Nasal/NP/ OP swab	DE ^[2]		DE ^[2]		14 July 2021
Jiangsu Diagnostics Biotechnology Co., Ltd	COVID-19 Antigen Rapid Test Cassette (Colloidal Gold)	1920	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 100%	97.58 % sensitivity 100 % specificity Nasal/NP/ OP swab	DE ^[2]		DE ^[2]		14 July 2021
Jiangsu Medomics medical technology Co.,Ltd.	SARS-CoV-2 antigen Test Kit (LFIA)	2006	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 94,1% at <Ct25) + Manufacturer specificity: 99,51%	Clinical Sensitivity: 97.73 % Clinical Specificity: 99.51 % Anterior nasal swab, NP swab	DE ^[2]		DE ^[2]		7 July 2021
Joinstar Biomedical	COVID-19 Rapid Antigen	1333	<i>Study design unknown</i>	96.1% sensitivity	AT, DE ^[2] ,		DE ^[2]		17 February

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance Data by manufacturer	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
Technology Co. Ltd	Test (Colloidal Gold)		DE: 96.1% sensitivity, 98.1% specificity SI: 96.1% sensitivity, 98.1% specificity, NP swab	98.1% specificity Nasal swab	PT, SI				2021
JOYSBIO (Tianjin) Biotechnology Co., Ltd.	SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold immunochromatography)	1764	<i>Prospective clinical field study</i>	98.13% sensitivity Nasal swab	AT, CZ, SI		CZ, DE ^[2] CH		10 May 2021
			FIND Evaluation CH (11 Feb 2021) 265 samples, Nasal swab Clinical sensitivities: - Days ≤ 7: 74.2%; - Ct ≤ 33: 78.9%; - Ct ≤ 25: 91.3%; Clinical specificity: 99.1%						
			<i>Study design unknown</i>						
			CZ: Meets the minimum performance requirements – see report for details.						
Labnovation Technologies Inc.	SARS-CoV-2 Antigen Rapid Test Kit	1266	<i>Study design unknown</i>	NP/OP swab	DE ^[2] , IT, SI		DE ^[2]		10 May 2021
			DE: 96.3% sensitivity, 97.3% specificity SI: 96.3% sensitivity, 97.3% specificity, NP/OP swab						
Lumigenex (Suzhou) Co., Ltd	PocRoc® SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)	2128	<i>Retrospective in vitro studies</i>	93.33% sensitivity 99.16% specificity Nasal/NP/OP swab	DE ^[2]		DE ^[2]		10 May 2021
			DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 99,16%						
LumiQuick Diagnostics Inc.	QuickProfile™ COVID-19 Antigen Test	1267	<i>Study design unknown</i>		BE, DE ^[2] , FR, SI,		DE ^[2]		10 May 2021
			BE: 94% sensitivity, 99% specificity, NP swab DE: 93.7% sensitivity, 98.8% specificity SI: 93.7% sensitivity, 98.8% specificity, NP swab						
LumiraDX	LumiraDx SARS-CoV-2 Ag Test	1268	<i>Study design unknown</i>	97.6% sensitivity 96.6% specificity Nasal swab	DE ^[2] , ES, SI	CH	DE ^[2] , ES, SKUP CH		17 February 2021
			DE: 93.8% sensitivity, 98.8% specificity SI: 97.6% sensitivity, 97.7% specificity, NP/Nasal swab SKUP/2021/124: 90% sensitivity, 97,8% specificity, NP swab (Scandinavian evaluation of laboratory equipment for point of care testing))						
MEDsan GmbH	MEDsan SARS-CoV-2	1180	<i>Study design unknown</i>	92.5% sensitivity	AT, BE, DE ^[2]	CH	DE ^[2]		17 February

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance Data by manufacturer	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
	Antigen Rapid Test		BE: 92.5% sensitivity, 99.8% specificity, Nasal/OP swab DE: 92.5% sensitivity, 99.8% specificity	99.8% specificity NP/OP swab			CH		2021
Merlin Biomedical (Xiamen) Co., Ltd.	SARS-CoV-2 Antigen Rapid Test Cassette	2029	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 90% at <Ct30 and 100% at <Ct25)	95.05% sensitivity 98.99% specificity Nasal/NP swab	DE ^[2]		DE ^[2]		16 June 2021
MEXACARE GmbH	MEXACARE COVID-19 Antigen Rapid Test	1775	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 99,1%	Clinical Sensitivity: 96.17 % Nasal swab	DE ^[2]		DE ^[2]		7 July 2021
möLab	mö-screen Corona Antigen Test	1190	<i>Study design unknown</i> DE: 97.25% sensitivity , 99.99% specificity	NP swab	DE ^[2] , IE		DE ^[2] , IE		10 May 2021
MP Biomedicals	Rapid SARS-CoV-2 Antigen Test Card	1481	<i>Study design unknown</i> BE: 96.4% sensitivity, 99% specificity, NP/OP swab DE: 96.39 % sensitivity, 99.03% specificity	96.17% sensitivity 99.16% specificity Nasal swab, Anterior nasal swab	AT, BE, DE ^[2]	CH	DE ^[2] CH		17 February 2021
Nal von minden GmbH	NADAL COVID -19 Ag +Influenza A/B Test	2104	<i>Study design unknown</i> DE: 97.6% sensitivity, 99.9% specificity	97% sensitivity 98% specificity NP swab	DE ^[2]		DE ^[2]		10 May 2021
Nal von minden GmbH	NADAL COVID -19 Ag Test	1162	<i>Prospective clinical field study</i> FIND evaluation CH (26 April 2021) 462 samples, NP swab Clinical sensitivities: - Days ≤ 7: 88.5%; - Ct ≤ 33: 92.4%; - Ct ≤ 25: 97.8%; Clinical specificity: 99.2% <i>Study design unknown</i> BE: 97.6% sensitivity, 99.9% specificity, NP/OP swab DE: 97.6% sensitivity, 99.9% specificity SI: 97.6% sensitivity, 99.9% specificity, NP/OP swab	97.6% sensitivity 99.9% specificity Nasal swab	AT, BE, CY DE ^[2] , FR, PT, SI		DE ^[2] , FR China	HR, SKUP	17 February 2021
NanoEntek	FREND COVID-19 Ag	1420	<i>Study design unknown</i>	94.12% sensitivity	DE ^[2]		DE ^[2]		10 May 2021

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance <i>Data by manufacturer</i>	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
			DE: 94.12% sensitivity , 100% specificity	100% specificity NP swab					
NanoRepro AG	NanoRepro SARS-CoV-2 Antigen Rapid Test	2200	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 94,1% at <Ct25) + Manufacturer specificity: 98.4%	97.2 % sensitivity 98.4% specificity Nasal/NP/OP swab	DE ^[2]		DE ^[2]		14 July 2021
NESAPOR EUROPA SL	MARESKIT COVID-19 ANTIGEN RAPID TEST KIT	Yes (2241)	<i>Study design unknown</i> ES: Independent validation study; Nasal test compared to nasal PCR. Sensitivity 95.24%, Specificity 100%.	Sensitivity: 95.24%, Specificity: 100% Nasal swab	ES		ES		23 July 2021
New Gene (Hangzhou) Bioengineering Co., Ltd.	COVID-19 Antigen Detection Kit	1501	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 92,5% at <Ct30 and 100% at <Ct25)	98% sensitivity Nasal swab	DE ^[2]		DE ^[2]		16 June 2021
Novatech	SARS-CoV-2 Antigen Rapid Test	1762	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 94,1% at <Ct25) + Manufacturer specificity: 100%	95 % sensitivity 100% specificity Nasal/ NP swab	DE ^[2]		DE ^[2]		14 July 2021
Oncosem Onkolojik Sistemler San. ve Tic. A.S.	CAT	1199	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 94,1% at <Ct25) + Manufacturer specificity: 98,04%	93.75% sensitivity 98.04% specificity Nasal swab	DE ^[2]		DE ^[2]		10 May 2021
PCL Inc.	PCL COVID19 Ag Rapid FIA	308	<i>Study design unknown</i> DE: 94,92 % sensitivity, 99,99 % specificity SI: 95.5% sensitivity, 98.6% specificity, NP/OP swab		FR, DE ^[2] , RO, SI		DE ^[2]		10 May 2021
PCL Inc.	PCL COVID19 Ag Gold	2243	<i>Prospective clinical field study</i> FR: Validation study data: 120 positive and 200 negative samples; sensitivity 92%, specificity: 100%		FR, PT		FR		7 July 2021
PerGrande Bio Tech Development Co., Ltd.	SARS-CoV-2 Antigen Detection Kit (Colloidal Gold Immunochromatographic Assay)	2116	<i>Study design unknown</i> DE: 94.28% sensitivity, 99.11% specificity	94.28% sensitivity 99.11% specificity NP/Nasal/OP swab	AT, DE ^[2]		DE ^[2]		10 May 2021
Precision Biosensor	Exdia COVI-19 Ag	1271	<i>Study design unknown</i>	93.9% sensitivity	SI, DE ^[2]	CH	DE ^[2]		17 February

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance <i>Data by manufacturer</i>	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
Inc.			DE: 93.88% sensitivity , 98% specificity SI:93.9% sensitivity, 98% specificity, NP swab	98% specificity NP swab			CH		2021
Prognosis Biotech	Rapid Test Ag 2019-nCov	1495	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 94,1% at <Ct25) + Manufacturer specificity: 99,58%	Clinical Sensitivity: 95.56 % Nasal swab	CY, DE ^[2]		DE ^[2]		7 July 2021
Qingdao Hightop Biotech Co. Ltd	SARS-CoV-2 Antigen Rapid Test (Immunochromatography)	1341	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct30 and 100% at <Ct25)	95% sensitivity Nasal swab	AT, DE ^[2]		DE ^[2]		17 February 2021
Quidel Corporation	Sofia SARS Antigen FIA	1097	<i>Study design unknown</i> BE: 96.7% sensitivity, 100% specificity, NP/nasal swab DE: 96.7% sensitivity , 100% specificity SI: 96.7% sensitivity, 100% specificity, NP/Nasal swab	96.7% sensitivity 100% specificity NP/Nasal swab	AT, BE, DE ^[2] , FI, NL ^[5] , PT, SI	CH	DE ^[2] , NL ^[5] CH	SI	17 February 2021
Rapid Pathogen Screening, Inc	LIAISON® Quick Detect Covid Ag Assay	Yes (2290)	<i>Study design unknown</i> IT: Independent validation study, 100 pos. and 100 neg. samples; sensitivity: 92.7% with Ct<25; specificity: 100%.	Sensitivity: 96.1%, Specificity: 97% NP and Nasal swab	IT		IT		23 July 2021
Roche (SD BIOSENSOR)	SARS-CoV-2 Rapid Antigen Test	1604	<i>Retrospective in vitro studies</i> FI: Validated in several laboratories (studies not published), meeting criteria. <i>Study design unknown</i> DE: 96.52% sensitivity, 99.68% specificity	96.52% sensitivity 99.2% specificity NP swab	AT, DE ^[2] , MT, NL, RO	CH, NO	DE ^[2] , FI		10 May 2021
Roche (SD BIOSENSOR)	SARS-CoV-2 Rapid Antigen Test Nasal	2228	<i>Prospective clinical field study</i> FIND evaluation DE (12 April 2021) 179 samples, nasal swab Clinical sensitivities: - Days ≤ 7: 81.2%; - Ct ≤ 33: 87.5%; - Ct ≤ 25: 100%; Clinical specificity: 99.3% Brazil (12 April 2021) 214 samples, nasal swab Clinical sensitivities: - Days ≤ 7: 81.2%;	Clinical Sensitivity: 89.6 % (Ct ≤ 30) 93.1 % (Ct ≤ 27) Clinical Specificity: 99.1 % Nasal swab	DK, SK	CH, UK	DE ^[2]		7 July 2021

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance <i>Data by manufacturer</i>	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
			- Ct ≤ 33: 91.7%; - Ct ≤ 25: 100%; Clinical specificity: 99.3% <i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 89.6% at <Ct30)						
Safecare Biotech (Hangzhou) Co. Ltd	COVID-19 Antigen Rapid Test Kit (Swab)	1489	<i>Study design unknown</i> DE: 97.27 % sensitivity , 99.42% specificity	97.04% sensitivity Nasal swab	AT, DE ^[2] , FR	CH	DE ^[2]		17 February 2021
Safecare Biotech (Hangzhou) Co. Ltd	Multi-Respiratory Virus Antigen Test Kit (Swab) (Influenza A+B/COVID-19)	1490	<i>Study design unknown</i> DE: 97.04% sensitivity , 99.44% specificity	97.04% sensitivity Nasal swab	DE ^[2]		DE ^[2]		10 May 2021
ScheBo Biotech AG	ScheBo SARS CoV-2 Quick Antigen	1201	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 95% at <Ct30 and 100% at <Ct25)	96.6% sensitivity (Ct ≤ 30) NP/ OP swab	DE ^[2]		DE ^[2]		16 June 2021
SD Biosensor Inc	STANDARD Q COVID-19 Ag Test Nasal	2052	<i>Prospective clinical field study</i> FIND evaluation DE (12 April 2021) 179 samples, nasal swab Clinical sensitivities: - Days ≤ 7: 81.2%; - Ct ≤ 33: 87.5%; - Ct ≤ 25: 100%; Clinical specificity: 99.3% Brazil (12 April 2021) 214 samples, nasal swab Clinical sensitivities: - Days ≤ 7: 81.2%; - Ct ≤ 33: 91.7%; - Ct ≤ 25: 100%; Clinical specificity: 99.3% <i>Retrospective in vitro studies</i>	Clinical Sensitivity: 97.12 % Clinical Specificity: 100 % Nasal swab	FI, PT, SK		DE ^[2] , FI, FR		7 July 2021

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance Data by manufacturer	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
			<p>FI: Validated in several laboratories (studies not published), meeting criteria.</p> <p><i>Study design unknown</i></p> <p>DE: Published study: https://www.medrxiv.org/content/10.1101/2021.01.06.20249009v1</p>						
SD BIOSENSOR Inc.	STANDARD F COVID-19 Ag FIA	344	<p><i>Prospective clinical field study</i></p> <p>FIND evaluation DE (10 Dec 2020) 676 samples, NP swab Clinical sensitivities: - Days \leq 7: 81.2%; - Ct \leq 33: 75%; - Ct \leq 25: 100%; Clinical specificity: 96.9%</p> <p>Brazil (10 Dec 2020) 453 samples, NP swab Clinical sensitivities: - Days \leq 7: 80.2%; - Ct \leq 33: 80.9%; - Ct \leq 25: 87.9%; Clinical specificity: 97.9%</p> <p>India (25 June 2020) 417 samples, NP swab Clinical sensitivities: - Days \leq 7: 61.8%; - Ct \leq 33: 53.6%; - Ct \leq 25: 68.5%; Clinical specificity: 99.5%</p> <p><i>Study design unknown</i></p> <p>BE: 96.5% sensitivity, 99.7% specificity, NP swab DE: 94% sensitivity 97% specificity</p>	94,09% sensitivity 98.52% specificity NP swab	AT, BE, BG, DE ^[2] , IT, LU, LV, NL ^[5] , PT, RO, SK	CH	DE ^[2] , IT, NL ^[5] , DK, CH, UK, BR	LU, PT	17 February 2021
SD BIOSENSOR Inc.	STANDARD Q COVID-19	345	<i>Prospective clinical field study</i>	96.52% sensitivity	AT, BE, BG,	ME, NO, CH	DE ^[2] , ES,	HR, IE, LU,	17 February

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance <i>Data by manufacturer</i>	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
	Ag Test		<p>FIND evaluation</p> <p>DE (10 Dec 2020) 1263 samples, NP swab Clinical sensitivities: - Days \leq 7: 80%; - Ct \leq 33: 87.8%; - Ct \leq 25: 100%; Clinical specificity: 99.3%</p> <p>Brazil (10 Dec 2020) 400 samples, NP swab Clinical sensitivities: - Days \leq 7: 90.7%; - Ct \leq 33: 91.9%; - Ct \leq 25: 95.9%; Clinical specificity: 97.6%</p> <p>CH (10 Dec 2020) 529 samples, NP swab Clinical sensitivities: - Days \leq 7: 89.8%; - Ct \leq 33: 91.8%; - Ct \leq 25: 97.2%; Clinical specificity: 99.7%</p> <p>India (22 April 2021) 334 samples, NP swab Clinical sensitivities: - Days \leq 7: 58.3%; - Ct \leq 33: 65.5%; - Ct \leq 25: 89.4%; Clinical specificity: 97.3%</p> <p>Peru (22 April 2021) 335 samples, NP swab Clinical sensitivities: - Days \leq 7: 81.4%; - Ct \leq 33: 83.3%; - Ct \leq 25: 96.2%; Clinical specificity: 99.6%</p> <p><i>Retrospective in vitro studies</i></p>	99.68% specificity NP swab	CY, DE ^[2] , DK, EE, ES, FI, FR, HR, IT, LU, LV, MT, NL ^[5] , PT, RO, SE, SK, SI		IT, NL ^[5] , DK, PT CH, UA, UK, BR, NO	SI, SE	2021

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance <i>Data by manufacturer</i>	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
			FI: Validated in several laboratories (studies not published), meeting criteria. <i>Study design unknown</i> BE: 96.5% sensitivity, 99.7% specificity, NP swab DE: 96.52% sensitivity, 99.68% specificity SI: 96.5% sensitivity, 99.7% specificity, NP swab						
SGA Medikal	V-Chek SARS-CoV-2 Ag Rapid Test Kit (Colloidal Gold)	1319	<i>Study design unknown</i> DE: 96.6% sensitivity, 99% specificity	96.6% sensitivity, Nasal swab	DE ^[2]		DE ^[2]		10 May 2021
SGA Medikal	V-Chek SARS-CoV-2 Rapid Ag Test (colloidal gold)	1357	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 94,1% at <Ct25) + Manufacturer specificity: 99,5%	Clinical Sensitivity: 96.60% Nasal swab	DE ^[2]		DE ^[2]		7 July 2021
Shenzhen Ultra-Diagnostics Biotech Co., Ltd	SARS-CoV-2 Antigen Test Kit	2017	<i>Study design unknown</i> BE: 92% sensitivity, 100% specificity, NP swab 100% sensitivity, 100% specificity, OP swab SI: 95.9% sensitivity, 99.9% specificity, NP/OP/Nasal swab	Clinical Sensitivity: 95.33 % (Nasal), 95.48(NP) Clinical Specificity: 99.16 % (Nasal), 99.61 % (NP)	AT, BE, ES, SI		BE, SI		10 May 2021
Shenzhen Lvshiyuan Biotechnology Co., Ltd.	Green Spring SARS-CoV-2 Antigen-Rapid test-Set	2109	<i>Study design unknown</i> DE: 98% sensitivity , 100% specificity	98% sensitivity 100% specificity NP/OP/Nasal swab	DE ^[2]		DE ^[2]		10 May 2021
Shenzhen Microprofit Biotech Co., Ltd	SARS-CoV-2 Antigen Test Kit (Colloidal Gold Chromatographic Immunoassay)	1967	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 100%	Clinical Sensitivity: 92.93 % Clinical Specificity: 100 % Nasal/NP/OP swab	DE ^[2] , ES		DE ^[2]		7 July 2021
Shenzhen Microprofit Biotech Co., Ltd.	SARS-CoV-2 Spike Protein Test Kit (Colloidal Gold Chromatographic Immunoassay)	1178	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 100%	Sensitivity: 86.3%, Specificity: 100% Nasal Swab	DE ^[2]		DE ^[2]		23 July 2021
Shenzhen Watmind Medical Co., Ltd	SARS-CoV-2 Ag Diagnostic Test Kit (Colloidal Gold)	1769	<i>Study design unknown</i> DE: 95.15% sensitivity , 99.12% specificity	95.15% Sensitivity Nasal swab	AT, DE ^[2] , FR		DE ^[2]		10 May 2021
Shenzhen Watmind	SARS-CoV-2 Ag Diagnostic	1768	<i>Retrospective in vitro studies</i>	Clinical Sensitivity:	DE ^[2]		DE ^[2]		7 July 2021

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance Data by manufacturer	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
Medical Co., Ltd	Test Kit (Immuno-fluorescence)		DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 99,13%	97.83 % (CT ≤ 33) Clinical Sensitivity: 90.08 % (Ct ≤ 36) Nasal swab					
Shenzhen Zhenrui Biotech Co., Ltd	Zhenrui ®COVID-19 Antigen Test Cassette	1574	<i>Study design unknown</i>						
			DE: 96% sensitivity 97% specificity	96% sensitivity Nasal swab	DE ^[2]		DE ^[2]		10 May 2021
Siemens Healthineers	CLINITEST Rapid COVID-19 Antigen Test	1218	<i>Study design unknown</i>						
			BE: 98.32% sensitivity, 99.6% specificity, NP swab 97.25% sensitivity, 100% specificity, Nasal swab SI: 96.7% sensitivity, 99.2% specificity, NP/Nasal swab	98.32% sensitivity (NP swab) 97.25% sensitivity 100% specificity (Nasal swab)	AT, BE, DE ^[2] , FR, HR, NL ^[5] , PT, SE, SI	CH	DE ^[2] , ES, NL ^[5]	HR, PT, SE ^[3]	17 February 2021
Sugentech, Inc.	SGTi-flex COVID-19 Ag	1114	<i>Retrospective in vitro studies</i>						
			DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct30 and 100% at <Ct25)	100% sensitivity 100% specificity OP/NP swab	AT, DE ^[2]		DE ^[2]		10 May 2021
TODA PHARMA	TODA CORONADIAG Ag	1466	<i>Study design unknown</i>						
			BE: 96.6% sensitivity, 100% specificity, NP/OP swab DE: 96.6% sensitivity, 100 specificity SI: 96.6% sensitivity, 100% specificity, NP/OP swab	98.6% sensitivity Nasal swab	BE, DE ^[2] , SI		DE ^[2]		10 May 2021
Triplex International Biosciences Co., Ltd	SARS-CoV-2 Antigen Rapid Test Kit	2074	<i>Retrospective in vitro studies</i>						
			DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 92,5% at <Ct30 and 100% at <Ct25)	98.33% sensitivity 100% specificity Nasal/OP/NP swab	DE ^[2]		DE ^[2]		16 June 2021
Triplex International Biosciences Co., Ltd, China	SARS-CoV-2 Antigen Rapid Test Kit	1465	<i>Retrospective in vitro studies</i>						
			DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 100%	98.51 % sensitivity Nasal swab	DE ^[2] , FR, PT		DE ^[2]		14 July 2021
Vitrosens Biotechnology Co., Ltd	RapidFor SARS-CoV-2 Rapid Ag Test	1443	<i>Retrospective in vitro studies</i>						
			DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct30 and 100% at <Ct25)						
			<i>Study design unknown</i>						
			SI: 97.3% sensitivity, 99% specificity, NP/OP/Nasal swab	97.3% sensitivity Nasal swab	DE ^[2] , SI		DE ^[2]		10 May 2021

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance Data by manufacturer	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
VivaChek Biotech (Hangzhou) Co., Ltd.	VivaDiag Pro SARS-CoV-2 Ag Rapid Test	2103	<i>Study design unknown</i> AT: 97,06% sensitivity, 100% specificity, all specimen types, i.e. N&OP&NP swab	97.04% sensitivity 99.9% specificity Nasal/OP/NP swab	AT, SI		AT, DE ^[2] , SI	AT	10 May 2021
Wuhan EasyDiagnosis Biomedicine Co., Ltd.	COVID-19 (SARS-CoV-2) Antigen-Test Kit	2098	<i>Study design unknown</i> DE: 96.15% sensitivity, 99.26% specificity	96.1% sensitivity 100% specificity Nasal/OP/NP swab	DE ^[2]		DE ^[2]		10 May 2021
Wuhan Life Origin Biotech Joint Stock Co., Ltd.	SARS-CoV-2 Antigen Assay Kit (Immunochromatography)	1773	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: xx%	92.67% sensitivity Nasal swab	DE ^[2]		DE ^[2]		14 July 2021
Wuhan UNscience Biotechnology Co., Ltd.	SARS-CoV-2 Antigen Rapid Test Kit	2090	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 99,57%	Sensitivity: 96.33% Specificity: 99.57% Nasal/NP/OP swab	DE ^[2]		DE ^[2] , FR		7 July 2021
Xiamen AmonMed Biotechnology Co., Ltd	COVID-19 Antigen Rapid Test Kit (Colloidal Gold)	1763	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 99.55%	93.2% sensitivity 99.55% specificity Nasal swab	DE ^[2]		DE ^[2]		10 May 2021
Xiamen Boson Biotech Co. Ltd	Rapid SARS-CoV-2 Antigen Test Card	1278	<i>Study design unknown</i> BE: 93.8% sensitivity, 100% specificity, NP swab DE: 96.49% sensitivity, 99.03% specificity	Not specified NP swab	AT, BE, BG, CY, DE ^[2] , FR, RO	CH	DE ^[2] CH		17 February 2021
Xiamen Wiz Biotech Co., Ltd	SARS-CoV-2 Antigen Rapid Test	1456	<i>Study design unknown</i> DE: 96.3% sensitivity, 100% specificity	96.3% sensitivity, Nasal swab	AT, DE ^[2]		DE ^[2]		10 May 2021
Xiamen Wiz Biotech Co., Ltd	SARS-CoV-2 Antigen Rapid Test (Colloidal Gold)	1884	<i>Study design unknown</i> DE: 95.91% sensitivity, 100% specificity	95.91% sensitivity 100% specificity Nasal swab	AT, DE ^[2]		DE ^[2]		10 May 2021
Zhejiang Anji Saianfu Biotech Co., Ltd	AndLucky COVID-19 Antigen Rapid Test	1296	<i>Study design unknown</i> DE: 97.5% sensitivity, 99.1% specificity	95.8% sensitivity, Nasal swab	AT, DE ^[2]		DE ^[2]		10 May 2021
Zhejiang Anji Saianfu Biotech Co., Ltd	reOpenTest COVID-19 Antigen Rapid Test	1295	<i>Retrospective in vitro studies</i> DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 94,1% at <Ct25) + Manufacturer specificity: 99%	95.8% sensitivity, Nasal swab	DE ^[2]		DE ^[2]		10 May 2021

Manufacturer	RAT commercial name	Device ID # ³	Clinical performance as reported by independent validation studies	Clinical performance <i>Data by manufacturer</i>	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	Included in EU common list since:
Zhejiang Orient Gene Biotech Co., Ltd	Coronavirus Ag Rapid Test Cassette (Swab)	1343	<i>Study design unknown</i>						
			BE: 98.32% sensitivity, 99.6% specificity, NP swab; 97.25% sensitivity, 100% specificity, Nasal swab DE: 96.72% sensitivity, 99.22% specificity	98.32 % sensitivity 99.6 % specificity Nasal/NP swab	AT, BE, BG, DE ^[2] , PT	CH, UK	DE ^[2]	SE ^[3]	17 February 2021
Zhuhai Lituo Biotechnology Co., Ltd.	COVID-19 Antigen Detection Kit (Colloidal Gold)	1957	<i>Retrospective in vitro studies</i>						
			DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 100%	96.12% sensitivity Nasal swab (CT _≤ 33); 99.59% sensitivity NP swab; 100% specificity Nasal swab (CT _≤ 33)	CZ, DE ^[2] , SI		DE ^[2]		14 July 2021